

# Projektplan

Tobias Nilsson

Version 1.0

Status:

Reviewed	Jesper Vesterberg	GRANSKATTDATUM
Approved	Krister Wiklund	GODKÄNTDATUM

Tobias Nilsson, toni0042  
Jesper Vesterberg, jeve0100

## Projektidentitet

Name	Ansvar	Telephone	e-mail
Tobias Nilsson	Documents		toni0042@student.umu.se
Jesper Vesterberg	Programmer		jeve0010@student.umu.se

**Client:** Institutionen för Fysik, Umeå universitet,  
Linnaeus väg 24,  
901 87  
Umeå.

**Contact:** Krister Wiklund, krister.wiklund@umu.se.

## Contents

0.1	Introduction . . . . .	1
<b>1</b>	<b>The set-up</b>	<b>1</b>
<b>2</b>	<b>Software</b>	<b>3</b>
<b>3</b>	<b>•</b>	<b>3</b>

## 0.1 Introduction

The laboratory *Electrical conductivity* in the courses *Solid State Physics* was upgraded from an old set-up using; a cryo system, with heat controller; a series of power supplies; a black box with measurement circuits; a bench multimeter and an ancient computer system running on MS-DOS. The upgrades consisted of removing the external power supplies, the black box and exchanging the computer for a newer laptop and the user interface, *Elledn*, was remade in Labview.

The purpose of the the laboratory is to investigate how the electrical conductivity of different materials behaves in the temperature region  $10 - 300$  K. The materials are a semi-conductor, InSb, a metal, Pt, and a super conductor,  $YBa_2Cu_3O_7$ .

## 1 The set-up

Equipment consist of:

- Cryogenic system:
  - Cooling chamber
  - Coarse vacuum pump
  - turbo molecular pump
  - Cryogenic system(compressor w. cooling system): HC-2, APD Cryogenics Inc.
- Controllers to cryogenic system:
  - Turbotronic NT10
  - Thermovac TM20
  - Scientific Instruments Inc. 9620-1 Silicon Diode
- Bench Multimeter: Keithley 2001 Multimeter
- Computer with Labview.

The materials are placed in the cooling chamber and the two vacuum pumps reduces the pressure to below  $10^{-5}$  Pascal. This reduces the amount of heat that needs to be removed to change the chambers temperature.

$\approx 10$  K is the lowest temperature for the cooling part of the cryogenic equipment, which removes heat at a constant pace. An external heater is

Table 1: *Pin layout for connection cables between cooling chamber and multimeter.*

25-pin D-sub	round 19-pin connector	Function			
1	A	Superconductor	Current	+	
2	B	"	"	-	
3	C	"	Voltage	+	
4	D	"	"	-	
5	E	Pt-100	Current	+	
6	F	"	"	-	
7	G	"	Voltage	+	
8	H	"	"	-	
9	J	InSb-plate	Current	+	
10	K	"	"	-	
11	L	"	Voltage	+	
12	M	"	"	-	
13	N	Temperature diode	Current	+	
14	P	"	"	-	
15	R	"	Voltage	+	
16	S	"	"	-	
17	T	free			
18	U	free			
25	V	Shield			

used to alter the temperature by adding heat to the system. This heater is controlled by *Scientific Instruments Inc. 9620-1 Silicon Diode*. Which in turn can be controlled by a computer through an GPIB-interface.

The controller determines the heating with the aid of an equation consisting of an integral and a differential term. To improve the performance and responsiveness of the cryogenic equipment in various ranges, these terms can be changed by the user. For more information see the *Scientific Instruments Inc. 9620-1 Silicon Diode manual*.

To determine the resistivity of the samples, each sample is connected to a 4 sense wire system. These sense wires are connected to a Keithley 2001 Multimeter through a switching card at the back. The multimeter is also equipped with a GPIB-interface, enabling computer controlled measurements.

The connections from the samples and an additional temperature diode, goes first through a round 19 pin contact to an 25 pin d-sub connector. See table 1 for pin layout on the cables.

## 2 Software

The software, ResSolidLabV1, is written in Labview and it can send the desired temperature to the

## 3 ●